

Instructions for Use



Revised version EN-BIO-UM-01.12.2016

Oosafe® IVF Laboratory Use Disinfectants

Ready-To-Use



Products

Oosafe® Disinfectant for CO₂ Incubators and Laminar Flow Hoods

Oosafe® Laboratory Surface Disinfectant

Oosafe® Hand Disinfectant

Oosafe® Disinfection Wipes

Attention!

- Check the expiration date on the label;
- Storage temperature is between 18 – 27°C;
- Acceptable temperature for transportation is from -18°C to 60°C;
- Avoid direct sunlight;
- Visually check for bottles to be undamaged and the sealing intact;
- Keep away from children;
- Check Material Data Sheet of the product for more information related to safety;
- Oosafe® disinfections bottles are made of high density polyethylene (HD-PE) and are recyclable. For the sake of the environment it is recommended to dispose bottles according to the local regulations.



- Cardboard boxes shall be disposed according to the local regulations for paper recycling;
- Plastic bottles and cardboard boxes are flammable. Keep away from fire.



First-aid measure!

- Inhalation: Supply fresh air, consult a doctor in case of complaints (applied only to Laboratory Surface Disinfectant).
- Skin contact: The product is not irritating for the skin.
- Eye contact: Rinse opened eye for 10 minutes under running water. Force to open eyelids. Consult a doctor if irritation persists.
- Ingestion: Rinse mouth and drink plenty of water. Do not give an unconscious person anything through the mouth. Do not cause vomiting. If symptoms persist, consult a doctor.

Additional information:

For more detailed information about the products, industrial safety and proper disposal refer to disinfectants safety data sheet, ingredient list and product specifications.

Oosafe® Disinfectant for CO₂ Incubators and Laminar Flow Hoods



Scope of application

Oosafe® Disinfectant for CO₂ incubators and laminar flow hoods is a transparent and odourless ready-to-use mixture. The disinfectant is alcohol-free, non-VOC releasing, non-fragranced, therefore it is a favourable disinfectant for medical equipment in IVF laboratories, including those with a low ventilation level. The disinfectant is effective against bacteria, yeasts, viruses, mycobactericidal (see table 1). The disinfectant does not cause corrosion of metals and does not affect organic soil.

Application

Recommended waiting time: 15 min.

Recommended use:

A method - Spray: Spray and distribute the disinfectant evenly on the surface with a sterile cloth. Wait 15 min for the disinfection effect. The surface is ready for work.

B method - Wipes: Wipe the surface evenly and wait 15 min for the disinfection effect.

Properties

- transparent and odourless;
- alcohol-free;
- ready-to-use;
- non-VOC releasing;
- effective against mycobacteria;
- effective against enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A (H5N1/H1N1), Ebola, also Adenovirus and Poliovirus;
- does not cause corrosion of metals;
- does not affect organic soil.

Intended use

- CO₂ incubators;
- Laminar flow hoods;
- Acrylic;
- metal;
- hi-Macs surfaces;
- glass;
- plastic.

Durability

After opening Oosafe® Disinfectant bottles or containers, the mixture maintains stability for 3 years under proper storage conditions.

Table 1: Efficacy table for Disinfectant for CO2 Incubators and Laminar Flow Hoods

Activity Spectrum	Method*	Contact Time
Bactericidal		
<i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i>	EN 13727+A1	5 min
	EN 1040	15 min
	EN 14561	15 min
Fungicidal		
<i>Candida albicans</i>	EN 1275	15 min
	EN 13624	15 min
	EN 14562	15 min
Mycobactericidal and Tuberculocidal		
<i>Mycobacterium avium</i> , <i>Mycobacterium terrae</i>	EN 14348	15 min
	EN 14563	15 min
Virucidal		
<i>Vaccinia virus</i> **	EN 14476+A1	15min
<i>Poliovirus</i>	EN 14476	15 min
<i>Adenovirus</i>	EN 14476+A1	5 min

*Efficacy is tested for clean conditions

**Effectiveness against all enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A (H5N1/H1N1), Ebola

Packaging Type

Packaging Type	Product Code
1. 50ml, spray bottle	OODIH-0050
2. 1L, spray bottle, individually packed in the carton box	OODIH-1000
3. 5L, refill container, individually packed in the cardboard protection	OODIH-5000
4. Wipes, bucket with 70 wipes	OODW-70

Oosafe® Laboratory Surface Disinfectant



Scope of application

Oosafe® Surface Disinfectant is a colourless, with a slightly smell of soap, the disinfectant meant to be diluted with water (1:50) before use. The disinfectant is alcohol-free, non-VOC releasing, non-fragranced, made for the disinfection and cleaning of laboratory surfaces and floors in IVF laboratories, including those with a low ventilation level. The non-ionic soap content helps cleaning and disinfecting at the same time. The disinfectant is effective against bacteria, fungi (candida) and viruses (see table 2).

Application

Recommended waiting time: 30 min.

Recommended use:

Method - Dilute Oosafe Surface Disinfectant 1:50 with water. Clean the surface and wait 30 min for the disinfection effect. If it is used for laboratory floor, please remember to rinse the floor 30 minutes after application, otherwise it may be slippery.

Properties

- clear-slightly yellow, and slightly aromatic;
- alcohol-free;
- non-VOC releasing;
- effective against bacteria, fungi (candida) and viruses;
- does not cause corrosion of metals;
- does not affect organic soil.

Intended use

- tiles;
- wood;
- plastic;
- rubber;
- glass;
- metal.

Durability

After opening Oosafe® Disinfectant bottles or containers, the mixture maintains stability for 3 years under proper storage conditions.

Table 2: Efficacy table for Laboratory Surface Disinfectant

Activity Spectrum	Method*	Contact Time
Bactericidal		
<i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i> , <i>Enterococcus hirae</i>	EN 14561	15 min
<i>Pseudomonas aeruginosa</i> , <i>Staphylococcus aureus</i>	EN 1040	15 min
Fungicidal		
<i>Candida albicans</i>	EN 1275 EN 14562 EN 13624	15 min
Virucidal		
<i>Vaccinia virus</i> *	EN 14476+A1	15 min

*Efficacy is tested for clean conditions

**Effectiveness against all enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A (H5N1/H1N1), Ebola

Packaging Type

Packaging Type	Product Code
1. 50 ml, bottle - sample	OODSF-00050
2. 2L, bottle with dose cup, individually packed in the carton box	OODSF-02000
3. 10L, refill container, individually packed in the cardboard protection	OODSF-10000

Oosafe® Hand Disinfectant



Scope of application

Oosafe® Hand Disinfectant is a colourless, odourless, ready to use water based disinfectant for hands. The disinfectant is alcohol-free, non-VOC releasing, non-fragranced, made for the use in IVF laboratories, including those with a low ventilation level.

Application

Recommended waiting time: 30 sec.

Recommended use:

Method – Rub hands with 3ml disinfectant for at least 30 sec. until hands are completely dry.

Properties

- clear, odourless;
- ready-to-use;
- alcohol-free;
- non-VOC releasing;
- effective against bactericidal, fungi (Candida);
- effective against enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A (H5N1/H1N1), Ebola;
- does not cause corrosion of metals;
- does not affect organic soil.

Intended use

- hands.

Intended use

After opening Oosafe® Disinfectant bottles or containers, the mixture maintains stability for 3 years under proper storage conditions.

Table 3: Efficacy table for Hand Disinfectant

Activity Spectrum	Method*	Contact Time
Bactericidal		
<i>Pseudomonas aeruginosa, Staphylococcus aureus</i>	EN 1040	30 secs
<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae</i>	EN 13727+A2 EN 14561	30 secs
<i>Escherichia coli</i>	EN 1500	30 secs
Fungicidal		
<i>Candida albicans</i>	EN 1275	30 secs
	EN 13624	30 secs
	EN 14562	30 secs
Virucidal		
<i>Vaccinia virus**</i>	EN 14476+A1	30 secs

*Efficacy is tested for clean conditions

**Effectiveness against all enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A (H5N1/H1N1), Ebola

Packaging Type

Packaging Type	Product Code
1. 50 ml, spray bottle - sample	OODH-0050
2. 05 L, bottle with pump: individually packed in the carton box	OODH-0500
3. 5 L, refill container, individually packed in the cardboard protection	OODH-5000

Manufacturer:

SparMED ApS

Ryttermarken 2, DK-3520 Farum, Denmark

E-m@il: info@sparmed.dk

Tel.: +45 39 40 25 03

Website: www.sparmed.dk