

HANDLING PIPETTE

PIPETTE FOR HANDLING EMBRYOS AND OOCYTES

The handling pipette is a single-use IVF laboratory micropipette used for transferring oocytes and embryos from one dish to another.

FEATURES


- IVF laboratory micropipette for dish-to-dish transfer of oocytes and embryos
- will not scratch dishes
- available in three lengths: 45, 55 and 94 mm
- available in two packaging variations: single packed and packed by ten pieces



ORIGINAL PICTURE

GYNO INFO HANDLING PIPETTE

Length	45, 55 or 94 mm
Tube OD	1.25 mm
Tube ID	1.00 mm
Sterile	sterilized by gamma irradiation
Packaging	250/shipping box (1/pouch) 1250/shipping box (10/pouch)
Quality control	MEA and LAL
Order number	#3145 (45 mm, 1/pouch) #3100 (55 mm, 1/pouch) #3125 (94 mm, 1/pouch) #3145 TP (45 mm, 10/pouch) #3100 TP (55 mm, 10/pouch) #3125 TP (94 mm, 10/pouch)



GMP Belgium Handling Pipette 09/2016