

# MIRI<sup>®</sup> GA

Gas and Temperature Validation Unit

*NEXT GENERATION Desktop CO<sub>2</sub> / O<sub>2</sub> & Temperature Validation Unit*

MIRI<sup>®</sup> GA is a tabletop device intended to make external incubator validation easier and safer. It is capable of monitoring the temperature (PT1000 connector) & gas concentration, flow and pressure. It can validate up to 6 chambers simultaneously 24 hours a day. The MIRI<sup>®</sup> GA can connect to any brand of incubator and is a perfect accessory to MIRI<sup>®</sup> TL and MIRI<sup>®</sup> Multiroom Incubators.



## Key Features



Constantly validate up to 6 x CO<sub>2</sub> / O<sub>2</sub> incubators



Controllable flow rate



Monitor up to 6 x PT1000 sensors



6 ports for sequential gas samples



Gas feedback returns sampled gas to incubator or exhaust



Audible and visual alarms



## As validation unit, MIRI<sup>®</sup> GA can provide

- Parameter data logging
- Adjustable alarm thresholds

Connected to a PC, the logging software will document parameters and optionally send e-mail alarms if set up.



Want to know more about this product? Scan the QR Code

**Esco Micro Pte. Ltd.**

21 Changi South Street 1 • Singapore 486 777 • Tel +65 6542 0833  
medical@escolifesciences.com • esco-medical.com

**Esco Global Offices:** Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam